

Wingmaster Elbow Positioning System

Part No. #ET1000

CLEANING AND STERILIZATION RECOMMENDATIONS

MATERIALS COMPATIBILITY

- ◁ Cleaning solutions must be compatible with; aluminum and acetal.
- ◁ Do not place in milk bath or lubricant.

POINT OF USE – O.R. WIPEDOWN

- ◁ Remove excess debris and gross soil.

STERILE PROCESSING

Disassemble all user movable components.

- ◁ All components can be disassembled without the use of tools.

Recommended washer disinfectant instructions:

Note: Do not place in milk bath or lubricant.

- ◁ Rinse in water to remove blood and saline solutions.
- ◁ Soak the components in an enzyme solution per enzyme manufacturer's directions.
- ◁ Put through washer/disinfectant.
- ◁ Dry all surfaces and wrap. (Optional sterilization tray or basket may be used.)

Recommended hand cleaning instructions:

Note: Do not place in milk bath or lubricant.

- ◁ Rinse in water to remove blood and saline solutions.
- ◁ Soak the components in an enzyme solution per enzyme manufacturer's directions.
- ◁ Wash and rinse the components.
- ◁ Inspect for remaining debris.
- ◁ Soak or wipe in low level disinfectant per hospital protocol.
- ◁ Rinse the components in warm water.
- ◁ Wash the components with an instrument detergent with a pH between 7 and 9.
- ◁ Rinse the components in warm water.
- ◁ Dry all surfaces and wrap. (Optional sterilization tray or basket may be used.)

STEAM AUTOCLAVE

Components may be double or two ply wrapped together or separately.

Cycle	Exposure Time 132°C (270°F)	Exposure Time 135°C (275°F)	Dry Time
Pre-Vacuum	4 Minutes		Follow wrap or container IFU
Gravity Displacement	20 Minutes	15 Minutes	Follow wrap or container IFU

Positioner must be allowed to cool to room temperature prior to use. Do not allow any component at an elevated temperature to contact the patient.