

Broadhead Pegboard Pegs

Lateral positioning pegs designed to reduce pressure.



- Pegs have a completely wipeable surface with no cracks or lamination.
- Broadhead pegs can be machine washed.
- Pegs are radiolucent to eliminate shadowing and make placement easier.
- More available peg length options than competitive products.
- Pegs and accessories compatible with Broadhead lateral, spine, and trauma pegboards.



4" PEG EXTENSIONS
#HP4020



6" PIVOT PEGS
#P6030



6" FIXED PEGS
#HP6020




6" PIVOT/FIXED COMBO
#HP6021



9" FIXED PEGS
#HP9020



9" PIVOT PEGS
#HP9030



9" PIVOT/FIXED COMBO
#HP9021



12" FIXED PEGS
#HP1220



www.matchgrademed.com | sales@matchgrademed.com | (920) 486-1955

