


# Ladderloc ZR Knee Support System

The familiar solution to total knee arthroplasty.



## Product Features

- Easy to use and easy to clean design.
- Familiar function with positive stop flexion adjustment.
- Reduces hands on the patient and perioperative resource time.
- High strength precision foot holder, will not bend, crack, or loose shape.
- Foot holder shape and construction allows for firm wrap and greater foot stabilization.
- Aircraft grade material for greater strength and durability.
- Respects the sterile field with Fieldloc technology.
- Longer and more stable rail system.

## System Includes


- Ladderloc Foot Holder (#LL916)
- Ladderloc ZR Rail (#LZR110)
- Universal Mounting Platform (#FL600)

## Dimensions

Adjustment Rail: 21”L x 7”W x 2.75”H

Foot Holder: 13”H x 10”L x 6.5”W

## ADDITIONAL PRODUCTS



**STERILE FOOT HOLDER PADS**  
#FP2500



**TRACTION TOWER SET**  
#LL2100



**PIVOT FOOT HOLDER**  
#LL920



[www.matchgrademed.com](http://www.matchgrademed.com) | [sales@matchgrademed.com](mailto:sales@matchgrademed.com)

