




Featured Product Showcase

2023

matchgrademed.com




Sterile Foot Holder Pads #FP2500

The patented Ladderloc Support System is an upgrade from a familiar design that is more versatile, easier to use and easier to clean.

This surgical knee positioner is made from a high strength aluminum alloy. It's longer, more stable and most importantly, respects the integrity of the sterile field.

Thanks to our patented Fieldloc Technology, the Ladderloc mounts securely on top of the surgical drapes and locks



into place for the procedure. This allows surgeons to have complete control of the flexion and extension of the patient's leg while providing an easy to use, solid base on which to operate.

The system reduces hands on the patient and preoperative resource time. Our foot holder is much stronger and won't lose its shape, bend, or crack even with heavy use. The hour glass shape and construction allows for a firm wrap that stays in place while offering greater foot stabilization.

Ladderloc System #LL2000



Ladderloc Rail
Ladderloc Standard Foot Holder
Fieldloc Mounting Platform




Fixed & Pivot Foot Holders #LL916 / #LL920



Fieldloc Universal Mounting Platform #FL600

Ladderloc Pivot #LL603



Ladderloc Sterile Component Set with Pivot Foot Holder

Solorail System

#SR2813

Solorail Rail
Solorail Pivot Foot Holder
Fieldloc Mounting Platform



The patented Solorail Support System is the premiere solution to total knee arthroplasty and the only full range surgical knee positioner that truly respects the integrity of the sterile field.

Our patented Fieldloc Technology mounts securely under the surgical drapes allowing the Solorail to be locked into place inside the sterile field. This allows surgeons to have complete control of the flexion and extension of their patient's leg while providing an easy to use, solid base on which to operate.

The smooth travel rail system provides stable and solid positioning with adjustability to optimize knee position for greater access and precision. The system reduces hands on the patient and preoperative resource time. The easy release adjustment block allows leg manipulation without changing the flexion adjustment, as well as reduces retractor stress on bones by allowing unrestricted anterior tibia translation.

matchgrademed.com

Solorail Deluxe

#SR4010

Solorail Rail
Solorail Pivot Foot Holder
Fieldloc Mounting Platform
Traction Tower Set




Fieldloc Universal Mounting Platform

#FL600

Pivot Foot Holder

#SR920




Sterile Foot Holder Pads

#FP2500

Case contains: 8 Sets of Foot Holder Pads + 3M Coban™


Made in the USA



Standard
#HP1000

Extended
#HP3000

Bariatric
#HP2000



The Broadhead Pegboard System is more than just a lateral hip positioning system. This unique patented system is the most versatile, easy to use, and easy to clean pegboard on the market.

The pegs are designed to eliminate pressure points, offer radiolucency without sacrificing strength and are machine washable.

Broadhead Pegboards have tapered ends and double wide pegs with nearly 4 inches of surface area. Both are designed specifically to better protect the

patient and reduce the possibility of pressure injuries occurring.

The pegs are fully radiolucent to allow for comprehensive patient imaging.

Our pegs are made of a premium strength polymer and come in multiple lengths and shapes. This provides more secure positioning of your patient as well as an array of positioning options during a procedure.

The Broadhead Mini is an ideal solution for lower extremity arthroscopic procedures.

Broadhead Pegboard System

#HP1000

Standard Broadhead Pegboard
9" Broadhead Pegs (4)
6" Broadhead Pegs (2)
4" Broadhead Peg Extensions (2)
Hold Down Brackets (2)

Also available as:
Bariatric Pegboard
#HP2000

Extended Leg Pegboard ABMS
#HP3000



6" Pegs
#HP6021

Sold in pairs:

Fixed - #HP6020
Pivot - #HP6030
Combo - #HP6021



9" Pegs
#HP9021

Sold in pairs:

Fixed - #HP9020
Pivot - #HP9030
Combo - #HP9021



12" Pegs
#HP1220

Sold in pairs:

Fixed - #HP1220



Mini Pegboard
#MPB1000

Mini Pegboard
7" Fixed Peg (1)
7" Pivot Peg (1)
Mounting Bracket

Made in the USA

matchgrademed.com


Speedloc Ankle Distractor

#ST1200

Speedloc Base
Speedloc Bracket
Speedloc Lever

The most innovative ankle distractor available. The Match Grade Medical Speedloc Ankle Distractor System is the ideal choice for arthroscopy and trauma. Precision and simplicity shine through in setup, operation, and cleaning.

Speedloc is quick and easy to adjust, giving the surgeon total control and feel of the traction being applied. The innovative, compact design eliminates the cables and knobs used in other outdated distractors, and brings you closer to your patient, with nothing to get in your way.



Rail Extension 18"

#EXT1018

Also sold as 24" - #EXT1024



Ankle Distractor Straps

#SS2000

Sterile ankle straps,
compatible with most ankle
distraction systems.

Case contains: 8 sterile
ankle straps

Wingmaster Deluxe Elbow Tower

#ET1000

The Wingmaster Elbow Tower merges precision and flexibility like never before. Our unique design gives you the ability to make adjustments when and where you need them. Ideal for distal humerus and olecranon fractures.

Includes standard Wingmaster system as well as extension support.

Uses standard round post rail clamp.



Adjustable Rolling Arm Table


#RC100



Wingmaster

Made in the USA


Speedloc



Lateral Positioning Pads

#FP4500

5 Pad Sets
Each set includes:
Body Pad (1)
Peg Pads (8)



Foot Holder Pads

#FP2500


Case contains: 8 sterile foot holder pads + 3M Coban™

Ankle Distractor Straps

#SS2000

Sterile ankle straps, compatible with most ankle distraction systems.

Case contains: 8 sterile ankle straps



Ladderloc Sterilization Tray

#LSC100

Universal Rail Clamp

#RC100

- Autoclavable
- Compatible with US & International Rails
- Accepts round or blade posts
- Can be completely disassembled for cleaning



Tibia/Femur Triangles

THREE SIZES

11" - TF011
14" - TF014
16" - TF016



Lateral Positioning System

#BLP2002

9" Fixed Positioning Peg (4)
Fixed Adjustment Base (4)
Standard Positioning Bracket (2)
Extended Positioning Bracket (2)
Offset Positioning Bracket (2)

Made in the USA

matchgrademed.com



Contact us for quotes, samples, and free trials.



matchgrademed.com | 920-486-1955