

Surgical Rail Clamp

Universal, Autoclavable Rail Clamp.



SURGICAL RAIL CLAMP #RC100

Product Features

- Fully autoclavable and designed for use in the sterile or non sterile field.
- Designed to work with both round and blade style brackets.
- Precision clamp provides secure connection over multiple layers of drape or directly on the table rail.
- Can be completely disassemble for cleaning.
- No hidden pivot points or springs to collect debris.
- A clamp truly designed to be used in the sterile field.

Specifications

- Fits North American and European Table Rails
- Accepts accessories with round posts up to 21/32" (18mm) diameter and flat bars up to 1" x 3/8" (26mm x 10mm).

ADDITIONAL PRODUCTS



**SPEEDLOC ANKLE
DISTRACTION SYSTEM**
#ST1200



STERILE ANKLE STRAP
#SS2000

PRODUCT IN USE



SURGICAL RAIL CLAMP
#RC100



www.matchgrademed.com | sales@matchgrademed.com

