


Broadhead Hip Positioning System

The complete solution to lateral hip positioning



STANDARD BROADHEAD PEGBOARD SYSTEM #HP1000

Product Features


- Broad contact surface for reducing the occurrence of pressure injuries.
- Pegs are Radiolucent allowing for comprehensive patient imaging.
- Tapered pegboard ends help reduce patient pressure points.
- Patent pending Double-Wide Pegs offer nearly 4 inches of flat bearing surface for unmatched support.
- Pegs are 100% wipeable, and machine washable.
- Multiple peg lengths and shapes provide an array of positioning options.
- Practical alternative for lower extremity arthroscopic thigh bolster and bumps.

System Includes


- 1 Standard Pegboard
- 9" Double-wide pegs (set of 4)
- 6" Double-wide pegs (set of 2)
- 4" Double-wide peg extensions (set of 2)
- Hold-down bracket (set of 2)

Disposable foam pad set also available for system.

ADDITIONAL PRODUCTS



12" BARIATRIC PEGS
#HP1220
(12" X 3.75")



BARIATRIC PEGBOARD
#HP2000
(47.5" X 22.75")



PEG EXTENSIONS
#HP4020
(4" X 3.75")



HOLD DOWN CLAMP SET
#HPC4000



www.matchgrademed.com | sales@matchgrademed.com

