

Solorail Deluxe Knee Support System



The complete solution to total knee arthroplasty



SOLORAIL DELUXE #SR4000

Product Features

- Easy to use and easy to clean design.
- Reduces hands on the patient and perioperative resource time.
- Applications for Total Joint, Sports Medicine, and MIS techniques.
- Precision adjustment through full flexion and extension.
- Wipeable surfaces without hidden crevices to hide organic debris.
- Lateral & medial support maintained during flexion adjustment.
- Reduces retractor stress on bones by allowing unrestricted anterior tibia translation.
- Thigh bolster gives additional lateral support.
- Foot holder pivots while under traction, improving access.


System Includes (#SR4000)

- Sterile Top Rail
- Universal Mounting Platform (#FL500)
- Pivot Foot Holder (#SR910)
- Adjustment Block
- Traction Tower Set (#SR2100)
 - Sterile Thigh Bolster
 - Distraction Post
 - Rail Base

INCLUDED PRODUCTS



UNIVERSAL MOUNTING PLATFORM #FL500



ADJUSTMENT BLOCK #SR280



PIVOT FOOT HOLDER #SR910



TRACTION TOWER SET #SR2100



www.matchgrademed.com | sales@matchgrademed.com

