


# Surgical Positioning Pads

The better padding alternative.



## Product Features

- Non-Sterile, disposable multi-layer deluxe lateral positioning foam pad.
- Individual peg pads for better patient support and reducing the occurrence of pressure injuries.
- Eliminates possibility for hidden organic debris that can be found on reusable pads.
- No concern about cleaning and storing heavy and easily damaged gel pads.
- Clean alternative for gel pads.

## System Includes

Each set includes 1 body cushion and 8 peg pads. Disposable Foam Pads available in a Case of 5 sets.

Case Includes:


- 5 - Body Cushions
- 40 - Peg Pads

## ADDITIONAL INFO



**STANDARD BODY PAD**  
(47" X 15" X 1")

- Top layer is a soft pliable urethane foam that adds gentle support and non-abrasive contact with the skin.
- The base layer is a lightly compressible supporting foam that will not bottom out and helps to eliminate impingement points.



**INDIVIDUAL PEG PADS**  
(3.5" X 6")

- These unique pads stick directly to hip positioning pegs with a temporary adhesive.
- The pads are easily placed on pegs before insertion into pegboard.
- The soft foam stays in place, will not wrinkle and will not bottom out.



[www.matchgrademed.com](http://www.matchgrademed.com) | [sales@matchgrademed.com](mailto:sales@matchgrademed.com)

