


Surgical Positioning Peg Pads

The better padding alternative.



Product Features

- Individual peg pads for better patient support and reducing the occurrence of pressure injuries.
- Eliminates possibility for hidden organic debris that can be found on reusable pads.
- These unique pads stick directly to hip positioning pegs with a temporary adhesive.
- The pads are easily placed on pegs before insertion into pegboard eliminating wrinkles that can concentrate pressure.
- The soft foam stays in place, will not wrinkle and will not bottom out.

System Includes

Each set includes 40 peg pads.
(5 sets of 8 pads)

ADDITIONAL PRODUCT



STANDARD BODY PAD (47" X 15" X 1")

- Top layer is a soft pliable urethane foam that adds gentle support and non-abrasive contact with the skin.
- The base layer is a lightly compressible supporting foam that will not bottom out and helps to eliminate impingement points.



www.matchgrademed.com | sales@matchgrademed.com

